

# SPE-Überziehschuh mit rutschfester Sohle

Einweg-Schuhüberzug aus Spunbond-Polypropylen, mit rutschfester Sohle. Zum Tragen über Ihren Schuhen in ungefährlichen Umgebungen zum Schutz vor Schmutz und Dreck. Ideal für unkritische Bereiche oder für Eingangsräume vor dem Einkleiden.

- SPP-Gewebe
- Länge gedehnt: 44 cm (A)
- Höhe: 16cm (B)
- Gewicht: 38 gsm
- Blau/weiß
- PSA-Kategorie I
- Nicht steril



Kat.-Nr.	Beschreibung	Menge/VE
16814471	SPE-Überziehschuh mit rutschfester Sohle	300

Artikel-Spezifikationen	Testmethode	Wert
Gewicht (gsm)	ISO 3801	38+0.2
Reißfestigkeit	MD ISO 13934-1	>55
(N/5 cm)	CD ISO 13934-1	>45
Dehnungsrate (%)	MD ISO 13934-1	>30
	CD ISO 13934-1	>30

Verpackung	
Verpackungsmenge	100 Schuhüberzüge pro Beutel 3 Beutel pro Karton
Kartonmaße	46 x 25 x 20cm

Austria: fishersci.at Belgium: fishersci.be Denmark: fishersci.dk  
Germany: fishersci.de Ireland: fishersci.ie Italy: fishersci.it  
Finland: fishersci.fi France: fishersci.fr Netherlands: fishersci.nl  
Norway: fishersci.no Portugal: fishersci.pt Spain: fishersci.es  
Sweden: fishersci.se Switzerland: fishersci.ch UK: fishersci.co.uk

© 2020 Thermo Fisher Scientific Inc. All rights reserved.  
Trademarks used are owned as indicated at fishersci.com/trademarks.

 **fisher scientific**  
part of Thermo Fisher Scientific